

# INSTALLATION OF SOFTWARE PTX C3



PLEASE SELECT **OK**



CLICK THE BUTTON TO INSTALL SOFTWARE PTX C3



CLICK CONTINUE



CONFIRM ALL WITH YES UNTIL SOFTWARE IS INSTALLED COMPLETELY.




# WINDOWS 7

For Windows 7 there are two specifications necessary – the pictures are in German language – but the settings are the same in English.

Click at the symbol of software PTX C3  with right mouse button.

Open SETTINGS.



Click the following settings under COMPATIBILITY

1. Windows XP - and confirm with **V**
2. Deactivation of... - confirm with **V**

